TrapLiner™ catheter

TrapLiner™ catheter

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The TrapLiner catheter is a rapid exchange guide extension catheter with a guidewire trapping balloon. It is offered in sizes compatible with 6F, 7F, and 8F guide catheters. The device has a proximal stainless steel pushrod covered on the distal portion with a semi-circular polymer. The trapping balloon is located proximal to the semi-circular polymer for trapping guidewires against the inner wall of a guide catheter. The pushrod is followed distally by a hydrophilic-coated full-round guide extension section.


The TrapLiner catheter has two radiopaque platinum-iridium marker bands on the guide extension segment. The proximal marker band is located near the collar. The distal marker band is located on the distal tip. The device also has a single radiopaque gold marker under the proximal end of the balloon. All three radiopaque markers enable visibility while using standard fluoroscopic methods. The device has two positioning marks located at 95cm (single mark) and 105cm (double mark) from the distal tip, which are not visible under fluoroscopy.

The TrapLiner catheter has been sterilized with ethylene oxide.

STERILE | EO

SPECIFICATIONS

Model / Size	5566 / 6F	5567 / 7F	5568 / 8F
Minimum guide catheter size	6F (0.070" / 1.78mm I.D.)	7F (0.078" / 1.98mm I.D.)	8F (0.088" / 2.24mm I.D.)
Pushrod outer diameter	0.020" / 0.51mm	0.025" / 0.64mm	
Guide extension inner diameter	0.056" / 1.42mm	0.062" / 1.57mm	0.071" / 1.80mm
Tip outer diameter (X)	1.70mm / 0.067" (5.1F)	1.90mm / 0.075" (5.7F)	2.16mm / 0.085" (6.5F)
Maximum guide catheter size	8F (0.091" / 2.31mm I.D.)		
Nominal pressure rating	12atm / 1216kPa		
Rated burst pressure	14atm / 1419kPa		

INDICATIONS

The TrapLiner catheter is intended for use in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, to facilitate placement of interventional devices, and to facilitate the exchange of an interventional device while maintaining the position of a guidewire within the vasculature.

CONTRAINDICATIONS

The TrapLiner catheter is contraindicated in vessels smaller than its tip outer diameter and in vessels in the neurovasculature.

WARNINGS

The TrapLiner catheter is provided sterile for single use only. Reuse of single-use devices creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Never advance the TrapLiner catheter into a vessel without a leading guidewire as vessel damage may result.

Never advance the TrapLiner catheter into a vessel smaller than its tip outer diameter. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the TrapLiner catheter, withdraw the catheter until the pressure returns to normal.

Due to the size and non-tapered tip of the TrapLiner catheter, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Never advance, withdraw, or rotate an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in separation of the catheter or guidewire tip, other device damage, or vessel injury.

Do not perform high pressure injections. Catheter damage may result if injection pressure exceeds 300psi.

Do not exceed rated burst pressure, as balloon burst may occur.

Always allow for bleed-back after deflating the TrapLiner balloon prior to continuing with the procedure. Failure to do so may result in air being introduced into the vasculature potentially resulting in an air embolism.

PRECAUTIONS

Do not use the catheter if the packaging is damaged. A damaged package could result in a breach of sterility or device damage.

Inspect the catheter prior to use for any damage. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

When coiling the catheter for flushing or periprocedural storage, avoid coiling the guide extension and balloon section, as catheter damage may result.

Flush the catheter lumen thoroughly with heparinized saline prior to use to prevent clot formation.

Exercise care in handling the catheter during a procedure to reduce the possibility of accidental breakage or kinking. Do not apply torque to the catheter during delivery, as catheter damage may result.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response, as catheter damage or vessel injury may occur.

Never advance the TrapLiner catheter more than 10cm beyond the tip of the guide catheter, as the TrapLiner catheter may become lodged in the guide catheter making it difficult to remove.

Do not withdraw an undeployed stent back into the TrapLiner catheter or advance the TrapLiner catheter over an undeployed stent while it remains in the body, as it may result in dislodging the stent from its shaft resulting in vessel occlusion and/or embolism. Instead, simultaneously pull both the TrapLiner catheter and undeployed stent back into the guide and out of the body.

ADVERSE EFFECTS

The following complications are generally associated with catheterization procedures and may occur when using the TrapLiner catheter:

- Cardiac arrest
- Embolism
- Infection
- Inflammatory response
- Slow-flow / occlusion
- Stent dislodgement
- Thrombosis
- Vessel dissection
- Vessel perforation

CLINICAL PROCEDURE

The TrapLiner catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- 1x Catheter

Other items required but not provided:

- Guide catheter with an inner diameter large enough to accommodate the specific model of TrapLiner catheter in use (See **Specifications**)
- Y-adaptor with hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter $\leq 0.014"$ / 0.36mm
- Inflation device
- 3-way stopcock
- 50:50 solution by volume of contrast medium and saline (for balloon inflation)
- Sterile heparinized saline (for flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the TrapLiner catheter packaging for damage.
2. Using sterile technique, remove the catheter dispenser coil from the packaging and transfer it to the sterile field.
3. Remove the TrapLiner catheter from the dispenser coil and inspect it for any damage.
4. Thoroughly flush the guide extension lumen of the catheter with heparinized saline solution.
5. Wet the distal section of the catheter with heparinized saline to activate the catheter's hydrophilic coating.
6. Using an inflation device and 3-way stopcock, perform a standard negative preparation of the balloon. Remove the inflation device, leaving the 3-way stopcock in the closed position attached to the proximal hub of the TrapLiner catheter.

DEPLOYMENT PROCEDURE

Deploy the TrapLiner catheter according to the following steps:

1. Secure the previously placed guidewire. Backload the distal tip of the TrapLiner catheter onto the proximal end of the guidewire and advance through the hemostasis valve and into the guide catheter.

Teleflex®

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis, MN 55369 USA

(866) 246-6990

(866) 804-9881

2.

Under fluoroscopy, advance the TrapLiner catheter beyond the distal tip of the guide catheter and into the desired location within the vessel.
3.

If performing an interventional procedure, backload the interventional device over the existing guidewire and advance the device through the guide catheter and TrapLiner catheter and into the desired vascular space.
NOTE: If a second wire is used during the intervention and encounters resistance within the guide catheter, pull the wire back several centimeters and slowly readvance.
4.

If trapping of the guidewire within the guide catheter for a device exchange is desired, withdraw the previously placed catheter over the existing 0.014" guidewire until the distal tip of the catheter is proximal to the radiopaque marker under the TrapLiner balloon. Confirm position using fluoroscopy. Inflate the TrapLiner balloon using an inflation device and slowly remove the catheter out of the guide catheter.
5.

If desired, a new catheter can be loaded onto the 0.014" guidewire and advanced until the catheter tip approaches the TrapLiner balloon.
6.

Deflate the balloon under fluoroscopy and allow for bleed-back from the hemostasis valve prior to continuing with the interventional procedure.
7.

After completing the procedure, remove the TrapLiner catheter prior to removing the guide catheter from the vessel.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Vascular Solutions LLC warrants that the TrapLiner catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special or consequential damages arising from the use of the TrapLiner catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

TrapLiner is a trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

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