SuperCross™ microcatheters

English/Instructions For Use.....	1
Česky/Návod k použití.....	2
Dansk/Brugsanvisning.....	4
Nederlands/Gebruiksaanwijzing.....	5
Eesti/Kasutusjuhised.....	7
Suomi/Käyttöohjeet.....	8
Français/Mode d'emploi.....	9
Deutsch/Gebrauchsanweisung.....	11
Ελληνικά/Οδηγίες χρήσης.....	13
Magyar/Használati utasítás.....	14
Italiano/Istruzioni per l'uso.....	16
Latviski/Lietošanas instrukcija.....	17
Lietuvių/Naudojimo instrukcijos.....	19
Norsk/Bruksanvisning.....	20
Polski/Instrukcja obsługi.....	22
Português/Instruções de utilização.....	23
Русский/Инструкции по применению.....	25
Español/Instrucciones de uso.....	26
Svenska/Bruksanvisning.....	28
Türkçe/Kullanım Talimatları.....	29
Slovenčina/Návod na použitie.....	31
Tiếng Việt/Hướng dẫn sử dụng.....	32
Hrvatski/Upute za uporabu.....	34
Srpski/Uputstvo za upotrebu.....	35
Română/Instrucțiuni de utilizare.....	36
Українська/Інструкції для застосування.....	38
ภาษาไทย/คำแนะนำสำหรับการใช้งาน.....	40

SuperCross™ Microcatheter

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The SuperCross microcatheter is a single lumen catheter designed for use in the coronary and/or peripheral vasculature. The SuperCross microcatheter provides guidewire support during interventional procedures and allows for the exchange of one distally located guidewire for another while maintaining access to distal vasculature. The SuperCross microcatheter has a hydrophilic coating on the distal 40cm of the straight tip catheters and the distal 72cm of the angled tip catheters to enhance deliverability to the target vasculature.

The SuperCross microcatheter has a single radiopaque marker band or coil located 0.035" / 0.89mm, 0.026" / 0.65mm, or 0.069" / 1.75mm from the distal tip, respectively. The catheter has white positioning marks located at 95cm (single mark) and 105cm (double mark) from the distal tip. The proximal end of the catheter incorporates a strain relief and a luer-lock adapter for flushing. SuperCross Angled Tip microcatheters have an angled tip straightener placed proximal to the angled section of the catheter. The angled tip straightener is intended to temporarily straighten the angled catheter tip to facilitate the backloading of a guidewire into the distal tip of an angled catheter.

The SuperCross microcatheter has been sterilized with ethylene oxide.



STERILE EO

Specifications – SuperCross Microcatheter Straight Tip

Model Number	5300	5301
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Marker band from distal tip, "B"	0.035" / 0.89mm	
Proximal shaft outer diameter	0.033" / 0.84mm (2.5F)	
Distal tip outer diameter	0.024" / 0.61mm (1.8F)	
Minimum guide catheter inner diameter	0.042" / 1.07mm	


Specifications – SuperCross FT Microcatheter Straight Tip

Model Number	5340	5341
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Marker band from distal tip, "B"	0.035" / 0.89mm	
Proximal shaft outer diameter	0.033" / 0.84mm (2.5F)	
Distal tip outer diameter	0.024" / 0.61mm (1.8F)	
Minimum guide catheter inner diameter	0.042" / 1.07mm	


Specifications – SuperCross Microcatheter 45° Angled Tip

Model Number	5302	5303
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	2.4mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	


Specifications – SuperCross Microcatheter 90° Angled Tip

Model Number	5304	5305
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	3.7mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	


Specifications – SuperCross Microcatheter 120° Angled Tip

Model Number	5306	5307
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	3.4mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	

Specifications – SuperCross Microcatheter 90° Extended Angled Tip

Model Number	5308	5309
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.069" / 1.75mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	5mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	

INDICATIONS

The SuperCross microcatheter is intended to be used in conjunction with steerable guidewires to access discrete regions of the coronary and/or peripheral vasculature. It may be used to facilitate placement and exchange of guidewires and other interventional devices and to subselectively infuse/deliver diagnostic and therapeutic agents.

CONTRAINDICATIONS

The SuperCross microcatheter is contraindicated for high pressure injections and for use in the cerebral vasculature.

WARNINGS

The SuperCross microcatheter is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Do not exceed the maximum injection pressure of 300 psi. Catheter damage or vessel injury could result.

Teleflex®

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis, MN 55369 USA

(866) 246-6990
(866) 804-9881

If an appropriate guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter *in vivo*. Catheter rupture and arterial injury could result. Identify and resolve the cause of the blockage or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

Do not advance the SuperCross microcatheter ahead of the guidewire, as vessel perforation may result.

PRECAUTIONS

Do not use the SuperCross microcatheter if the packaging is damaged. A damaged package could indicate a breach of sterility or device damage.

Inspect the catheter prior to use. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

Care should be taken not to crush the catheter. Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the wire lumen and difficulty while exchanging wires.

The catheter lumen should be flushed with sterile, heparinized saline prior to use to ensure catheter is free from debris which could be introduced into the body resulting in an embolism.

Exercise care while handling the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response as catheter damage or vessel injury may occur.

Do not insert the SuperCross Angled Tip into the guide catheter without removing the angled tip straightener as it could prevent the hemostasis valve from sealing around the catheter.

Check that all fittings and attachments are secure and observe that the syringe is free of air bubbles, prior to injection, so that air is not introduced into the catheter during use. Failure to do so may result in air being introduced into the vasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the use of the SuperCross catheter include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel spasm

CLINICAL PROCEDURE

The SuperCross microcatheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- Single-use disposable SuperCross microcatheter
- Angled tip straightener for SuperCross Angled Tip microcatheters only
- Dispenser coil

Note: The SuperCross microcatheter straight tip dispenser coil includes a flushing luer.

Other items required but not provided:

- If using the SuperCross microcatheter straight tip, a guiding catheter with a minimum I.D. of 0.042" / 1.07mm fitted with a hemostasis valve (Tuohy-Borst type)
- If using the SuperCross microcatheter angled tip, a guiding catheter with a minimum I.D. of 0.056" / 1.42mm fitted with a hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter $\leq 0.014"$ / 0.36mm
- 10ml syringe for flushing
- Sterile heparinized saline

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the SuperCross microcatheter packaging and components for damage.
2. Utilizing sterile technique, remove the SuperCross microcatheter from its packaging and transfer the dispenser coil to the sterile field.

3. If using the SuperCross microcatheter straight tip, attach a 10ml syringe filled with sterile, heparinized saline to the flushing luer on the dispenser coil. Completely flush the coil to activate the catheter's hydrophilic coating.
4. If using the SuperCross microcatheter angled tip, remove the catheter from the dispenser coil and wet the distal section with heparinized saline to activate the catheter's hydrophilic coating. Advance the angled tip straightener towards the distal tip such that it straightens the tip angle and maintains the tip in this straightened position in preparation for guidewire loading.
5. Inspect the SuperCross microcatheter for any bends or kinks. Thoroughly flush the SuperCross microcatheter with heparinized saline.

DEPLOYMENT PROCEDURE

Deploy the SuperCross microcatheter according to the following steps:

1. Backload the distal tip of the SuperCross microcatheter onto the proximal end of the exchange length guidewire that is already in place in the distal vasculature.
2. If the angled tip straightener has been used to straighten the tip of an angled SuperCross microcatheter tip to facilitate guidewire backloading, ensure that the guidewire is backloaded past the proximal edge of the angled tip straightener
3. If the angled tip straightener has been used and once the guidewire has been backloaded as described in step 2, carefully tear the angled tip straightener off the distal end of the catheter by gripping the proximal tab on the angled tip straightener and peeling it in the direction of the distal tip of the microcatheter. Discard the angled tip straightener once it has been removed from the catheter.
4. Under fluoroscopy, carefully advance the catheter until the marker at the distal tip is visible in the desired vascular space.
WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter or vessel perforation.
5. If guidewire exchange is desired, after fluoroscopically confirming that the SuperCross microcatheter tip is positioned in the desired vascular space, slowly remove the guidewire.
6. Insert the desired exchange length guidewire into the luer-lock guidewire entry port. Advance the guidewire until it exits the catheter tip.
7. Under fluoroscopy, ensure that the guidewire is sufficiently advanced into the distal vasculature. Fix the position of the guidewire and slowly advance the SuperCross microcatheter into the desired vascular space.
8. If subselective infusion/deliver of diagnostic or therapeutic agents is desired, after fluoroscopically confirming that the SuperCross microcatheter tip is positioned in the desired vascular space, slowly remove the guidewire.
9. Create an air-free, fluid-to-fluid connection between the catheter hub and the injection syringe. Aspirate prior to injection, and observe the syringe to remove any air bubbles prior to injecting into the catheter.
10. Using the injection syringe, manually infuse the diagnostic or therapeutic agent.

WARNING: Do not exceed the maximum recommended injection pressure of 300 psi when injecting through the SuperCross catheter. Catheter damage or vessel injury could result.

11. Ensure proper positioning of the exchange length guidewire prior to removing the SuperCross catheter and proceeding with the procedure.
12. Dispose of SuperCross microcatheter following standard hospital procedure.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Vascular Solutions LLC warrants that the SuperCross microcatheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special, or consequential damages arising from the use of the SuperCross microcatheter. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of Vascular Solutions LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

SuperCross is a trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

See the International Symbols Glossary on page 43.

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