Twin-Pass™

dual access catheters

English/Instructions For Use.....	1
Česky/Návod k použití.....	2
Dansk/Brugsanvisning.....	3
Nederlands/Gebruiksaanwijzing.....	3
Eesti/Kasutusjuhised.....	4
Suomi/Käyttöohjeet.....	5
Français/Mode d'emploi.....	6
Deutsch/Gebrauchsanweisung.....	7
Ελληνικά/Οδηγίες χρήσης.....	8
Magyar/Használati utasítás.....	9
Italiano/Istruzioni per l'uso.....	10
Latviski/Lietošanas instrukcija.....	11
Lietuvių/Naudojimo instrukcijos.....	12
Norsk/Bruksanvisning.....	13
Polski/Instrukcja stosowania.....	13
Português/Instruções de utilização.....	14
Español/Instrucciones de uso.....	15
Svenska/Bruksanvisning.....	16
Türkçe/Kullanım Talimatları.....	17
Slovenčina/Návod na použitie.....	18
Hrvatski/Upute za uporabu.....	19
Srpski/Uputstvo za upotrebu.....	20
Română/Instrucțiuni de utilizare.....	20
Русский/Инструкция по применению.....	21
Українська/Інструкції для застосування.....	22

Teleflex®

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis, MN 55369 USA

(866) 246-6990
(866) 804-9881

Twin-Pass™ Dual Access Catheter

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Twin-Pass dual access catheter is a dual lumen catheter designed for use in the coronary and/or peripheral vasculature. The Twin-Pass catheter consists of an over-the-wire (OTW) lumen that runs the length of the catheter and a rapid exchange (RX) delivery lumen on the distal segment. The Twin-Pass catheter has a working length of 135cm and is compatible with various guidewire diameters (see specification table). The Twin-Pass catheter also has hydrophilic coating on the distal 18cm of the catheter.

The Twin-Pass catheter has white positioning marks located at 95cm (single mark) and 105cm (double marks) from the distal tip, respectively. The Twin-Pass catheter has a marker band located 1mm from the distal tip and a second marker band located 11mm from the distal tip to identify the distal end of the OTW lumen.


The Twin-Pass catheter has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Model Number	5200
Max. guidewire diameter - RX lumen	0.014" / 0.36mm
RX lumen inner diameter	0.017" / 0.43mm
Max. guidewire diameter - OTW lumen	0.014" / 0.36mm
OTW lumen inner diameter	0.018" / 0.46mm
Proximal shaft outer diameter	0.038" / 0.97mm (2.9F)
Distal shaft outer diameter	0.047" / 1.19mm (3.5F)
Distal tip outer diameter	0.026" / 0.66mm (2F)
Minimum guide catheter inner diameter	0.066" / 1.68mm

Model 5200

INDICATIONS

The Twin-Pass catheter is intended to be used in conjunction with steerable guidewires in order to access discrete regions of the coronary and peripheral arterial vasculature, to facilitate placement and exchange of guidewires and other interventional devices, and for use during two guidewire procedures.

The Twin-Pass catheter is also used to sub selectively infuse/deliver diagnostic or therapeutic agents.

CONTRAINDICATIONS

The Twin-Pass catheter is contraindicated for high pressure injections and for use in the cerebral vasculature.

WARNINGS

The Twin-Pass catheter is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Do not advance the Twin-Pass catheter without a guidewire in place through an RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in vessel injury.

Do not exceed the maximum recommended infusion pressure of 300 psi (2068 kPa) when injecting through the Twin-Pass catheter as catheter damage or vessel injury could result.

If a guidewire of appropriate diameter cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter damage or vessel injury could result. Identify and resolve the cause of the blockage, or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in catheter damage or vessel injury.

PRECAUTIONS

Do not use the Twin-Pass catheter if the packaging has been damaged. A damaged package could indicate a breach of sterility or device damage.

Inspect the catheter prior to use. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumens must be flushed with sterile, heparinized saline prior to use to ensure catheter is free from debris which could be introduced into the body.

Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the guidewire lumen, difficulty while inserting the catheter or guidewires, or pressure damping.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response as catheter damage or vessel injury may occur.

Check that all fittings and attachments are secure and observe that the syringe is free of air bubbles, prior to injection, so that air is not introduced into the catheter during use. Failure to do so may result in air being introduced into the vasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the use of the Twin-Pass catheter include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel spasm

CLINICAL PROCEDURE

The Twin-Pass catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- Single-use Twin-Pass catheter

Other items required but not provided:

- Guide catheter with an inner diameter large enough to accommodate the specific model of Twin-Pass catheter in use
- Appropriate guidewires
- 10ml syringe (for system flushing)
- Sterile heparinized saline (for system flushing)
- Transducer (if measuring pressure)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the Twin-Pass catheter packaging and components for damage.
2. Using sterile technique, transfer the dispenser coil with the Twin-Pass catheter into the sterile field.
3. Attach a 10ml syringe filled with sterile, heparinized saline to the luer-lock guidewire entry port of the Twin-Pass catheter and thoroughly flush the catheter.
4. Flush the Twin-Pass catheter with sterile, heparinized saline to activate the hydrophilic coating.
5. Remove the Twin-Pass catheter from the dispenser coil and inspect it for any bends or kinks.

DEPLOYMENT PROCEDURE

Deploy the Twin-Pass catheter according to the following steps:

1. Backload the rapid exchange segment of the Twin-Pass catheter onto the guidewire that is already in place in the distal vasculature.
WARNING: Do not advance the Twin-Pass catheter without a guidewire in place through an RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in vessel injury.
2. Under fluoroscopy, carefully advance the Twin-Pass catheter into the desired distal vascular space.
WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in catheter damage or vessel injury.
3. To deliver an exchange-length guidewire through the OTW lumen:
 - Load the guidewire into the proximal lumen of the Twin-Pass catheter.
 - Advance the guidewire until it exits the OTW lumen into the distal vascular space.
 - Fix both guidewires using standard guidewire exchange techniques and carefully withdraw the Twin-Pass catheter until

- the distal tip exits the hemostatic valve and both wires can be secured.
- To deliver fluid through the OTW lumen:
 - Create an air-free, fluid-to-fluid connection between the OTW lumen and the injection syringe.
 - Aspirate prior to injection, and observe the syringe to remove any air bubbles prior to injecting into the catheter.
 - Using the injection syringe, manually infuse the diagnostic or therapeutic agent.

WARNING: Do not exceed the maximum recommended infusion pressure of 300 psi (2068 kPa) when injecting through the Twin-Pass catheter. Catheter damage or vessel injury could result.

Precaution: Check that all fittings and attachments are secure so that air is not introduced into the catheter during the procedure.

- Dispose of Twin-Pass catheter following standard hospital procedure.

STORAGE AND HANDLING

No special or handling conditions.

LIMITED WARRANTY

Vascular Solutions LLC warrants that the Twin-Pass catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special or consequential damages arising from the use of the Twin-Pass catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

Twin-Pass is a trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

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