Twin-Pass™ TORQUE

dual access catheter

English/Instructions for Use.....	1
Česky/Návod k použití.....	2
Dansk/Brugsanvisning.....	2
Nederlands/Gebruiksaanwijzing.....	3
Eesti/Kasutusjuhised.....	4
Français/Mode d'emploi.....	5
Deutsch/Gebrauchsanweisung.....	6
Ελληνικά/Οδηγίες χρήσης.....	7
Magyar/Használati utasítás.....	7
Italiano/Istruzioni per l'uso.....	8
Latviski/Lietošanas instrukcija.....	9
Lietuvių/Naudojimo instrukcijos.....	10
Norsk/Bruksanvisning.....	11
Polski/Instrukcja użytkowania.....	11
Português/Instruções de utilização.....	12
Español/Instrucciones de uso.....	13
Svenska/Bruksanvisning.....	14
Türkçe/Kullanım Talimatları.....	15
Slovenčina/Návod na použitie.....	15
Suomi/Käyttöohjeet.....	16
Hrvatski/Upute za uporabu.....	17
Srpski/Uputstvo za upotrebu.....	18
Română/Instrucțiuni de utilizare.....	19
Русский/Инструкция по применению.....	19
Українська/Інструкції для застосування.....	20

Twin-Pass Torque dual access catheter

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Twin-Pass Torque dual access catheter ("Twin-Pass") is a dual lumen catheter consisting of an over-the-wire (OTW) lumen that runs the full length of the catheter and a rapid exchange (RX) delivery lumen on the distal segment. The distal section of the catheter has a hydrophilic coating.



The catheter has two radiopaque platinum-iridium marker bands to identify the two exit ports. The device has two positioning marks located at 95cm (single mark) and 105cm (double marks) from the distal tip, which are not visible under fluoroscopy.

There are no diagnostic or therapeutic agents known to be incompatible with the Twin-Pass catheter. The Twin-Pass catheter has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Model Number	5201
Working length	135cm
Compatible guidewire diameter	≤ 0.014" / 0.36mm
OTW lumen exit port deflection angle	10°
Compatible guide catheter size	≥ 5F (≥ 0.056" / 1.42mm ID)
Dual-lumen outer diameter	1.17mm / 0.046" (3.5F)
Distal tip outer diameter	0.71mm / 0.028" (2.1F)
Maximum allowable outer diameter (dual-lumen)	1.30mm / 0.051" (3.9F)

INDICATIONS

The Twin-Pass catheter is intended to access discrete regions of the coronary and/or peripheral vasculature. It may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.

CONTRAINDICATIONS

None

WARNINGS

The Twin-Pass catheter is provided sterile for single use only. Reuse of single-use devices creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

If an appropriately-sized diameter guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter rupture and arterial injury may result. Identify and resolve the cause of the blockage or replace the catheter with a new one.

Never advance, withdraw, or rotate an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in separation of the catheter or guidewire tip, other device damage, or vessel injury.

Do not advance the catheter without an in-place guidewire through the RX lumen. Advancement of the catheter without a leading guidewire may result in vessel injury.

Do not perform high pressure injections through the catheter. Catheter damage may result if injection pressure exceeds 300 psi.

PRECAUTIONS

Do not use the Twin-Pass catheter if the packaging is damaged. A damaged package could result in a breach of sterility or device damage.

Inspect the catheter prior to use for any kinks. Do not use a damaged catheter. An inability to advance or withdraw the catheter may occur.

Exercise care when handling the catheter during a procedure to reduce the possibility of accidental breakage, bending, or kinking.

Flush both catheter lumens thoroughly with heparinized saline prior to use to prevent clot formation.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing

the resultant tip response, as catheter damage or vessel injury may occur.

Care should be taken not to crush the catheter. Excessive tightening of the hemostatic valve onto the catheter shaft may result in damage to the guidewire lumen and difficulty when inserting the catheter or guidewires.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system.

This device has not been evaluated for use in the neurovasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the Twin-Pass catheter include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel Spasm

CLINICAL PROCEDURE

The Twin-Pass catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- 1x Catheter

Other items required but not provided:

- Guide catheter or sheath with an inner diameter large enough to accommodate the Twin-Pass catheter (See **Specifications**)
- Hemostasis valve (tuohy borst type)
- Guidewire(s) with diameter ≤ 0.014" / 0.36mm
- 10mL syringe (for system flushing)
- Sterile heparinized saline (for system flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the Twin-Pass catheter packaging for damage.
2. Using sterile technique, remove the catheter from the packaging and transfer it to the sterile field.
3. Use a 10mL syringe filled with sterile, heparinized saline to thoroughly flush both catheter lumens.
4. Wet the outside of the catheter with sterile, heparinized saline to activate the hydrophilic coating.

DEPLOYMENT PROCEDURE

Deploy the Twin-Pass catheter according to the following steps:

1. Backload the distal tip of the Twin-Pass catheter onto the proximal end of the already positioned guidewire and advance through the hemostasis valve.
2. Carefully advance the Twin-Pass catheter to the desired location in the vasculature.
3. To deliver a second guidewire through the OTW lumen:
 - Load the guidewire into the proximal hub of the catheter and advance under fluoroscopy until the distal tip of the guidewire reaches the proximal radiopaque marker on the catheter.
 - Slowly advance the guidewire out of the OTW lumen exit port and identify the direction of wire deflection using fluoroscopy.
 - Rotate the catheter as needed in order to direct the guidewire tip towards the desired location for advancement. Advance the guidewire to the desired location in the vasculature.
4. To deliver fluid through the OTW lumen:
 - Create an air-free, fluid-to-fluid connection between the catheter hub and the injection syringe.
 - Aspirate prior to injection to remove any air bubbles in the catheter.
 - Using the injection syringe, manually deliver the diagnostic or therapeutic agent.
5. In order to remove the catheter, fix the guidewire(s) using standard catheter exchange techniques and carefully withdraw the catheter until the distal tip exits the hemostatic valve and the guidewire(s) can be secured.
6. Dispose of catheter following standard hospital procedure.

STORAGE & HANDLING

No special storage or handling conditions.

Teleflex®

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis, MN 55369 USA

(866) 246-6990

(866) 804-9881

LIMITED WARRANTY

Vascular Solutions LLC warrants that the Twin-Pass dual access catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special or consequential damages arising from the use of the Twin-Pass dual access catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intl

Twin-Pass is a trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

See the International Symbols Glossary on page 23.

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