R350™ guidewire

English/Instructions For Use.....	1
Dansk/Brugsanvisning.....	1
Nederlands/Gebruiksaanwijzing.....	2
Eesti/Kasutusjuhised.....	2
Français/Mode d'emploi.....	3
Deutsch/Gebrauchsanweisung.....	4
Ελληνικά/Οδηγίες χρήσης.....	4
Italiano/Istruzioni per l'uso.....	5
Latviski/Lietošanas instrukcija.....	5
Lietuvių/Naudojimo instrukcijos.....	6
Norsk/Bruksanvisning.....	6
Polski/Instrukcja użytkowania.....	7
Português (Brasil)/Instruções de uso.....	7
Português (Portugal)/Instruções de utilização.....	8
Español/Instrucciones de uso.....	9
Svenska/Bruksanvisning.....	9
Türkçe/Kullanım Talimatları.....	10
Slovenčina/Návod na použitie.....	10
Suomi/Käyttöohjeet.....	11
Čeština/Návod k použití.....	11
Magyar/Használati utasítás.....	12
Hrvatski/Upute za uporabu.....	12
Srpski/Uputstvo za upotrebu.....	13
Română/Instrucțiunile de utilizare.....	13
Русский/Инструкция по применению.....	14
Українська/Інструкції для застосування.....	15

R350™ guidewire

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.


DEVICE DESCRIPTION

The R350 guidewire is a 350cm nitinol mandrel guidewire with a 0.013" (0.33mm) diameter and a straight distal tip. It has a 5cm gold-plated tungsten coil on the distal end that is visible under fluoroscopic methods. The proximal 150cm has a PTFE coating, and the distal 200cm has a hydrophilic coating.

The R350 guidewire has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

INDICATIONS

The R350 guidewire is indicated for use in percutaneous procedures to introduce and position catheters and other interventional devices within the coronary and/or peripheral vasculature.

CONTRAINDICATIONS

The R350 guidewire is contraindicated for use in the neurovasculature.

WARNINGS

The R350 guidewire is provided sterile for single use only. Reuse of single-use devices creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the guidewire against resistance may result in separation of the guidewire tip, damage to the guidewire, or vessel perforation.

PRECAUTIONS

Do not use the R350 guidewire if the packaging is damaged. A damaged package could result in a breach of sterility or device damage.

Inspect the guidewire prior to use for any kinks. Do not use a damaged guidewire. Vessel damage and/or inability to advance or withdraw the guidewire may occur.

Exercise care in handling the guidewire during a procedure to reduce the possibility of accidental breakage, bending, or kinking.

When the R350 guidewire is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the guidewire without observing the resultant tip response, as guidewire damage or vessel injury may occur.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the R350 guidewire include, but are not limited to, the following:

- cardiac tamponade
- embolism
- infection
- hematoma
- inflammatory response
- thrombosis
- vessel dissection
- vessel perforation
- vessel spasm

CLINICAL PROCEDURE

The R350 guidewire should be used by physicians trained on the procedures for which it is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- 1x R350 guidewire

Other items required but not provided:

- Guide catheter
- Hemostasis valve (Tuohy-Borst type)
- Sterile heparinized saline (for flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the R350 guidewire packaging and components for damage.
2. Using sterile technique, remove the guidewire from its packaging and transfer it to the sterile field.

3. Thoroughly flush the guidewire with heparinized saline solution.

DEPLOYMENT PROCEDURE

Deploy the R350 guidewire according to the following steps:

1. Following placement of a guide catheter into the vasculature, insert the distal tip of the guidewire into the hemostasis valve attached to the guide catheter. Deliver the guidewire through the hemostasis valve and into the guide catheter.
2. Under fluoroscopy, advance the guidewire out of the guide catheter and into the selected vessel.
3. If required, secure the guidewire and backload a device over the guidewire and into the vasculature.
4. Once the intervention is complete, maintain guidewire position while withdrawing the device over the guidewire prior to removing the guidewire from the vasculature.
5. Dispose of the R350 guidewire following standard hospital procedures.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Teleflex Medical LLC warrants that the R350 guidewire is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Teleflex Medical LLC to be defective in workmanship or materials. Teleflex Medical LLC shall not be liable for any incidental, special or consequential damages arising from the use of the R350 guidewire. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Teleflex Medical LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Teleflex Medical LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF TELEFLEX MEDICAL LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

R350 is a trademark of Teleflex Innovations S.à.r.l., or Teleflex Medical, each a part of Teleflex Incorporated.

See the International Symbols Glossary on page 18.

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