WARRIOR™ guidewire

English / Instructions For Use.....	1
Español / Instrucciones de uso.....	2
Português (Brasil) / Instruções de uso.....	2
Nederlands / Gebruiksaanwijzing.....	3
Français / Notice d'utilisation.....	4
Deutsch / Gebrauchsanweisung.....	5
Italiano / Istruzioni per l'uso.....	5
Polski / Instrukcja użycia.....	6
Tiếng Việt / Hướng dẫn sử dụng.....	7
Čeština / Návod k použití.....	8

Warrior™ guidewire

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Warrior guidewire is a 0.014" diameter stainless steel core guidewire with a 0.009" diameter tapered distal tip. It is available in 200cm and 300cm lengths. The distal 20cm of the guidewire has a spring coil, of which the distal 2.5cm is platinum and visible under fluoroscopic methods. The guidewire has a straight shapeable tip with a tip load of 14 grams. The distal 20cm of the guidewire has a hydrophilic coating and the proximal portion has a PTFE coating.


Model 7388 has a modified proximal end to allow for guidewire extension.

The Warrior guidewire has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Model	7388	7389
Length (x)	200cm	300cm
Shaft outer diameter	0.014"	
Tip outer diameter	0.009"	
Tip load	14.0g	

INDICATIONS

The Warrior guidewire is intended for use in percutaneous procedures to introduce and position catheters and other interventional devices within the coronary and peripheral vasculature.

CONTRAINDICATIONS

The Warrior guidewire is contraindicated for use in the neurovasculature.

WARNINGS

The Warrior guidewire is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the guidewire against resistance may result in separation of the guidewire tip, damage to the guidewire, or vessel injury.

The Warrior guidewire has a stiff distal tip which may present a higher risk of perforation or injury than a guidewire with a more flexible distal tip. Therefore, use a guidewire with the least stiff distal tip that can be delivered across the lesion, and use care to minimize the risk of perforation or other damage to vessels.

PRECAUTIONS

Do not use the guidewire if the packaging is damaged. A damaged package could indicate a breach of sterility or device damage.

Inspect the guidewire prior to use for any damage. Do not use a damaged guidewire. Vessel damage and/or inability to advance or withdraw the guidewire may occur.

Exercise care in handling the guidewire during a procedure to reduce the possibility of accidental breakage or kinking.

When the guidewire is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the guidewire without observing the resultant tip response, as guidewire damage or vessel injury may occur.

When shaping the guidewire tip, do not use an object with a sharp edge as this may damage the guidewire.

Always secure the guidewire prior to advancing or withdrawing an interventional device over it. Failure to do so may result in loss of guidewire position and/or vessel injury.

Before inserting the guidewire into an interventional device, always flush the guidewire with saline to ensure the hydrophilic coating is thoroughly wet. Failure to do so may inhibit wire movement within the interventional device.

For 200cm length only: Do not attach a torque device to the proximal 2cm of the guidewire shaft, as it may damage the extendable portion of the guidewire, preventing the ability to attach to the DOC® Guide Wire Extension device.

ADVERSE EFFECTS

Potential adverse effects (in alphabetical order) that may be associated with the Warrior guidewire include, but are not limited to, the following:

- Cardiac tamponade
- Embolism
- Hematoma
- Infection
- Inflammatory response
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel spasm

CLINICAL PROCEDURE

The Warrior guidewire should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- 1x Warrior guidewire

Other items required but not provided:

- Guide catheter
- Guidewire introducer
- Hemostasis valve
- Sterile heparinized saline (for flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the guidewire packaging and components for damage.
2. Utilizing sterile technique, remove the guidewire hoop from its packaging and transfer to the sterile field.
3. Thoroughly flush the guidewire contained within the hoop with heparinized saline to activate the hydrophilic coating.
4. If desired, the guidewire tip may be carefully shaped using standard techniques

PRECAUTION: When shaping the guidewire tip, do not use an object with a sharp edge as this may damage the guidewire.

DEPLOYMENT PROCEDURE

Deploy the Warrior guidewire according to the following steps:

1. Following placement of a guide catheter and/or other interventional catheter into the vasculature, insert the distal tip of the guidewire into the catheter hub or hemostasis valve attached to the catheter. A guidewire introducer will be required if delivering the guidewire into a hemostasis valve.
2. Advance the guidewire through the catheter and into the vasculature.
3. If desired, attach a torque device to the proximal end of the guidewire at any time.
4. If desired, catheter exchange over the guidewire can be performed. Secure the guidewire prior to removing any catheters over the guidewire or backloading any new catheters onto the guidewire.
5. Once the intervention is complete, remove the guidewire from the vasculature.
6. Dispose of the guidewire following standard hospital procedures.

Guidewire Extension (200cm model only):

If desired, the Abbott Vascular DOC Guide Wire Extension device may be utilized to facilitate catheter exchange over a standard length (200cm) guidewire. Follow the manufacturer's instructions in order to perform a catheter exchange using the extension device.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Teleflex Medical LLC warrants that the Warrior guidewire is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Teleflex Medical LLC to be defective in workmanship or materials. Teleflex Medical LLC shall not be liable for any incidental, special or consequential damages arising from the use of the Warrior guidewire. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Teleflex Medical LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Teleflex Medical LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF TELEFLEX MEDICAL LLC.

Teleflex™

Teleflex Medical LLC
3015 Carrington Mill Blvd
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

Warrior is a trademark of Teleflex Innovations S.à.r.l., or Teleflex Medical, each a part of Teleflex Incorporated.

DOC is a registered trademark of Abbott Cardiovascular Systems, Inc.

See the International Symbols Glossary on page 11

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DOC je registrovaná ochranná známka společnosti Abbott Cardiovascular Systems, Inc.