Turnpike™ catheters

Turnpike™ catheters

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

This device has not been evaluated for use in the neurovasculature.

DEVICE DESCRIPTION

The Turnpike catheters are single lumen catheters designed for use in the coronary and/or peripheral vasculature. The catheters' shafts are constructed of two polymer layers that encapsulate a braid and a dual-layer coil. The distal 60cm of each catheter has a hydrophilic coating.

The Turnpike is available in four versions: the Turnpike catheter (*long, flexible, radiopaque distal tip*), the Turnpike LP catheter (*lower profile & increased flexibility*), the Turnpike Spiral catheter (*adds an outer coil onto the distal shaft for additional rotational advancement*), and the Turnpike Gold catheter (*outer coil on the distal shaft and a gold-plated, stainless steel, threaded distal tip*).

There are no diagnostic or therapeutic agents known to be incompatible with the Turnpike catheters.

The Turnpike catheters have been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Figure 1: Turnpike Catheter (standard)


Table 1: Turnpike Catheter (standard)

Model number	5642	5643
Working length (X)	135cm	150cm
Compatible guidewire diameter	≤ 0.014" / 0.36mm	
Compatible guide catheter size	≥ 5F (≥ 0.056" / 1.42mm ID)	
Proximal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal shaft outer diameter	0.86mm / 0.034" (2.6F)	
Distal tip outer diameter	0.53mm / 0.021" (1.6F)	
Maximum allowable outer diameter	1.02mm / 0.040" (3.1F)	

Figure 2: Turnpike LP Catheter


Table 2: Turnpike LP Catheter

Model Number	5638	5639
Working length (X)	135cm	150cm
Compatible guidewire diameter	≤ 0.014" / 0.36mm	
Compatible guide catheter size	≥ 5F (≥ 0.056" / 1.42mm ID)	
Proximal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal shaft outer diameter	0.74mm / 0.029" (2.2F)	
Distal tip outer diameter	0.53mm / 0.021" (1.6F)	
Maximum allowable outer diameter	1.02mm / 0.040" (3.1F)	

Figure 3: Turnpike Spiral Catheter


Table 3: Turnpike Spiral Catheter

Model Number	5640	5641
Working length (X)	135cm	150cm
Compatible guidewire diameter	≤ 0.014" / 0.36mm	
Compatible guide catheter size	≥ 5F (≥ 0.056" / 1.42mm ID)	
Proximal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal tip outer diameter	0.53mm / 0.021" (1.6F)	
Maximum allowable outer diameter	1.02mm / 0.040" (3.1F)	

Figure 4: Turnpike Gold Catheter


Table 4: Turnpike Gold Catheter

Model Number	5621	5622
Working length (X)	135cm	150cm
Compatible guidewire diameter	≤ 0.014" / 0.36mm	
Compatible guide catheter size	≥ 5F (≥ 0.056" / 1.42mm ID)	
Proximal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal shaft outer diameter	0.97mm / 0.038" (2.9F)	
Distal tip outer diameter	0.71mm / 0.028" (2.1F)	
Maximum allowable outer diameter	1.07mm / 0.042" (3.2F)	

INDICATIONS

The Turnpike catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.

CONTRAINDICATIONS

The Turnpike Spiral and Turnpike Gold catheters are contraindicated for use in vessels with an effective diameter smaller than 1mm.

WARNINGS

The Turnpike catheters are provided sterile for single use only. Reuse of single-use devices creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

If an appropriately-sized guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter rupture and arterial injury may result. Identify and resolve the cause of the blockage or replace the catheter with a new one.

Never advance, withdraw, or rotate an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, other device damage, or vessel injury.

Always advance the Turnpike catheters over an existing in-place guidewire. Do not use the Turnpike catheters to cross a lesion or occlusion ahead of the distal tip of the guidewire, as vessel injury may result.

Do not perform high pressure injections. Catheter damage may result if injection pressure exceeds 300psi.

Do not rotate the catheter more than two (2) consecutive 360° rotations in either direction if the distal tip becomes lodged and is not also rotating, as it may result in separation of the catheter, damage

en	English/Instructions for Use	1
cs	Česky/Návod k použití	2
da	Dansk/Brugsanvisning	3
nl	Nederlands/Gebruiksaanwijzing	5
et	Eesti/Kasutusjuhised	6
fi	Suomi/Käyttöohjeet	7
fr	Français/Mode d'emploi	9
de	Deutsch/Gebrauchsanweisung	10
el	Ελληνικά/Οδηγίες χρήσης	12
hu	Magyar/Használati utasítás	13
it	Italiano/Istruzioni per l'uso	14
lv	Latviski/Lietošanas instrukcija	16
lt	Lietuvių/Naudojimo instrukcijos	17
no	Norsk/Bruksanvisning	18
pl	Polski/Instrukcja stosowania	19
pt	Português/Instruções de utilização	21
es	Español/Instrucciones de uso	22
sv	Svenska/Bruksanvisning	23
tr	Türkçe/Kullanım Talimatları	25
sk	Slovenčina/Návod na použitie	26
vi	Tiếng Việt/Hướng dẫn sử dụng	27
hr	Hrvatski/Upute za uporabu	29
sr	Srpski/Uputstvo za upotrebu	30
ro	Română/Instrucțiuni de utilizare	31
ru	Русский/Инструкция по применению	32
uk	Українська/Інструкція для застосування	34
th	ภาษาไทย/คำแนะนำสำหรับการใช้งาน	35

Teleflex™

Teleflex Medical LLC
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

to the catheter, or vessel injury. If using the Turnpike or Turnpike LP, unwind the catheter to its relaxed position and apply backward tension in order to dislodge the tip. If using the Turnpike Spiral or Turnpike Gold, rotate the catheter in the counter-clockwise direction with backward tension in order to dislodge the tip.

Never advance the Turnpike Spiral or Turnpike Gold catheters into a vessel with an effective diameter smaller than 1mm. The catheter could become lodged and result in vessel injury.

PRECAUTIONS

Do not use the Turnpike catheters if the packaging is damaged. A damaged package could result in a breach of sterility or device damage.

Inspect the catheter prior to use. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumen should be flushed with sterile, heparinized saline prior to use to ensure catheter is free from debris which could be introduced into the body resulting in an embolism.

Care should be taken not to crush the catheter. Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the wire lumen and/or difficulty while exchanging wires.

Exercise care in handling the catheter during a procedure to reduce the possibility of accidental breakage, bending, or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response, as catheter damage or vessel injury may occur.

If using the Turnpike Spiral or Turnpike Gold, the catheter may only be rotated in the counter-clockwise direction to remove. Vessel damage and/or inability to withdraw the catheter may occur.

Check that all fittings and attachments are secure and observe that the syringe is free of air bubbles, prior to injection, so that air is not introduced into the catheter during use. Failure to do so may result in air being introduced into the vasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the use of the Turnpike catheters include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Septal hematoma
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel spasm

CLINICAL PROCEDURE

The Turnpike catheters should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- 1x Catheter
- 1x Dispenser coil with flushing luer

Other items required but not provided:

- Guide catheter or sheath with an inner diameter large enough to accommodate the Turnpike catheter in use (See **Specifications** section)
- Hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter $\leq 0.014'' / 0.36\text{mm}$
- 10mL sterile syringe (for flushing the dispenser coil and the catheter lumen)
- Sterile heparinized saline (for system flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the catheter packaging for damage.
2. Utilizing sterile technique, remove the catheter dispenser coil from the packaging and transfer it to the sterile field.
3. Attach a 10mL syringe filled with sterile, heparinized saline to the luer-lock guidewire entry port of the catheter and thoroughly flush the lumen.
4. Attach a 10mL syringe filled with sterile, heparinized saline to the flushing luer of the dispenser coil and completely flush the coil.
5. Remove the catheter from the dispenser coil and inspect it for any bends or kinks.

DEPLOYMENT PROCEDURE

Deploy the Turnpike catheters according to the following steps:

1. Secure the previously placed guidewire. Backload the distal tip of the Turnpike catheter onto the proximal end of the guidewire and advance through the hemostasis valve.
2. Carefully advance the Turnpike catheter through the vasculature under fluoroscopy. If resistance is encountered while advancing the catheter through the vessel, rotate the catheter in either a clockwise or counter-clockwise direction while gently pushing on the catheter to advance. If using the Turnpike Spiral or Turnpike Gold, rotate the catheter in the clockwise direction only to advance.
3. If guidewire exchange is desired, after fluoroscopically confirming that the Turnpike catheter distal tip is positioned in the desired vascular space, slowly remove the guidewire and insert the new guidewire into the luer-lock guidewire entry port. Advance the guidewire until it exits the catheter tip. Under fluoroscopy, ensure that the guidewire is sufficiently advanced into the distal vasculature.
4. If subselective infusion/delivery of diagnostic or therapeutic agents is desired, after fluoroscopically confirming that the Turnpike tip is positioned in the desired vasculature space, slowly remove the guidewire.
5. Create an air-free, fluid-to-fluid connection between the catheter hub and the injection syringe. Aspirate prior to injection, and observe the syringe to remove any air bubbles prior to injecting into the catheter.
6. Using the injection syringe, manually infuse the diagnostic or therapeutic agent.
7. Fix the guidewire using standard catheter exchange techniques and carefully withdraw the Turnpike catheter. If resistance is encountered while removing the Turnpike catheter, rotate the catheter in either a clockwise or counter-clockwise direction while gently pulling on the catheter to remove. If using Turnpike Spiral or Turnpike Gold, rotate the catheter in the counter-clockwise direction only to remove.
8. Dispose of catheter following standard hospital procedure.

WARNING: Do not rotate the catheter more than two (2) consecutive 360° rotations in either direction if the distal tip is not also rotating and advancing, as it may result in separation of the catheter, damage to the catheter, or vessel injury.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Teleflex Medical LLC warrants that the Turnpike catheters are free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Teleflex Medical LLC to be defective in workmanship or materials. Teleflex Medical LLC shall not be liable for any incidental, special or consequential damages arising from the use of the Turnpike catheters. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Teleflex Medical LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Teleflex Medical LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF TELEFLEX MEDICAL LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

Turnpike is a registered trademark of Teleflex Innovations S.à.r.l., or Teleflex Medical LLC, each a part of Teleflex Incorporated.

See the International Symbols Glossary on page 37.

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