SuperCross™ microcatheters

SuperCross™ Microcatheters

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION AND PERFORMANCE CHARACTERISTICS

The SuperCross Angled Tip microcatheters are single lumen catheters designed for use in the coronary and/or peripheral vasculature. The SuperCross microcatheters provide guidewire support during interventional procedures and allows for the exchange of one distally located guidewire for another while maintaining access to distal vasculature. The SuperCross Angled Tip microcatheters have a hydrophilic coating on the distal 72cm to enhance deliverability to the target vasculature.

The SuperCross Angled Tip microcatheters have a single radiopaque coil located 0.026" / 0.65mm (standard tip) or 0.069" / 1.75 mm (extended tip) from the distal tip. The catheter has white positioning marks located at 95 cm (single mark) and 105 cm (double mark) from the distal tip. The proximal end of the catheter incorporates a strain relief and a luer-lock adapter for flushing. The SuperCross Angled Tip microcatheters have an angled tip straightener placed proximal to the angled section of the catheter. The angled tip straightener is intended to temporarily straighten the angled catheter tip to facilitate the backloading of a guidewire into the distal tip of an angled catheter.


The SuperCross Angled Tip microcatheters have been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS


Specifications – SuperCross Microcatheters 45° Angled Tip

Model Number	5302	5303
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	2.4mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	

Specifications – SuperCross Microcatheters 90° Angled Tip


Model Number	5304	5305
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	3.7mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	

Specifications – SuperCross Microcatheters 120° Angled Tip


Model Number	5306	5307
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.026" / 0.65mm	

Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)
Distal tip outer diameter	0.031" / 0.79mm (2.4F)
Angled tip length	3.4mm
Minimum guide catheter inner diameter	0.056" / 1.42mm

Specifications – SuperCross Microcatheters 90° Extended Angled Tip

Model Number	5308	5309
Max. guidewire diameter	0.014" / 0.36mm	
Working length, "A"	130cm	150cm
Radiopaque coil length, "B"	0.400" / 10mm	
Radiopaque coil from distal tip, "C"	0.069" / 1.75mm	
Proximal shaft outer diameter	0.042" / 1.07mm (3.2F)	
Distal tip outer diameter	0.031" / 0.79mm (2.4F)	
Angled tip length	5mm	
Minimum guide catheter inner diameter	0.056" / 1.42mm	

INTENDED PURPOSE

The SuperCross Angled Tip microcatheters are intended for use in percutaneous procedures to facilitate the placement and exchange of interventional devices and to deliver diagnostic and therapeutic agents.

INDICATIONS FOR USE

The SuperCross microcatheter is intended to be used in conjunction with steerable guidewires to access discrete regions of the coronary and/or peripheral vasculature. It may be used to facilitate placement and exchange of guidewires and other interventional devices and to subselectively infuse/deliver diagnostic and therapeutic agents.

CONTRAINDICATIONS

The SuperCross microcatheter is contraindicated for high pressure injections and for use in the cerebral vasculature.

PATIENT TARGET GROUPS

The SuperCross microcatheters are intended to be used in adult patients undergoing percutaneous catheterization procedures.

INTENDED USERS

The SuperCross microcatheter should be used by physicians trained on the procedures for which the device is intended (i.e. Interventional Cardiologists). The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

CLINICAL BENEFITS TO BE EXPECTED

The ability to facilitate access to discrete regions of the coronary and peripheral vasculature through a single puncture site for the placement and exchange of interventional devices and the delivery of diagnostic and therapeutic agents.

WARNINGS

The SuperCross microcatheter is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Do not exceed the maximum injection pressure of 300 psi. Catheter damage or vessel injury could result.

If an appropriate guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter rupture and arterial injury could result. Identify and resolve the cause of the blockage or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

en	Instructions For Use	1
bg	Инструкции за употреба	2
hr	Upute za uporabu	4
cs	Návod k použití	5
da	Brugsanvisning	7
nl	Gebruiksaanwijzing	8
fi	Käyttöohjeet	9
fr	Notice d'utilisation	11
de	Gebrauchsanweisung	12
el	Οδηγίες χρήσης	14
hu	Használati utasítás	15
it	Istruzioni per l'uso	17
lt	Naudojimo instrukcija	18
no	Bruksanvisning	20
pl	Instrukcja użycia	21
pt	Instruções de utilização	22
ro	Instrucțiuni de utilizare	24
ru	Инструкция по применению	25
sr	Uputstvo za upotrebu	27
sk	Návod na použitie	28
sl	Navodila za uporabo	30
es	Instrucciones de uso	31
sv	Bruksanvisning	32
th	คำแนะนำสำหรับการใช้งาน	34
tr	Kullanma Talimatı	35
uk	Інструкції для застосування	36
vi	Hướng dẫn sử dụng	38

Teleflex Medical LLC
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA
USA: (866) 246 6990
International: +1 919 544 8000

EU Authorized
Representative and Importer:

EC REP

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

CH REP

Swiss AR Services AG
Industriestrasse 47
CH-6300 Zug

Teleflex™
www.teleflex.com

Do not advance the SuperCross microcatheter ahead of the guidewire, as vessel perforation may result.

PRECAUTIONS

Do not use the SuperCross microcatheter if the packaging is damaged. A damaged package could indicate a breach of sterility or device damage.

Inspect the catheter prior to use. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

Care should be taken not to crush the catheter. Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the wire lumen and difficulty while exchanging wires.

The catheter lumen should be flushed with sterile, heparinized saline prior to use to ensure catheter is free from debris which could be introduced into the body resulting in an embolism.

Exercise care while handling the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response as catheter damage or vessel injury may occur.

Do not insert the SuperCross Angled Tip into the guide catheter without removing the angled tip straightener as it could prevent the hemostasis valve from sealing around the catheter.

Check that all fittings and attachments are secure and observe that the syringe is free of air bubbles, prior to injection, so that air is not introduced into the catheter during use. Failure to do so may result in air being introduced into the vasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the use of the SuperCross catheter include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel spasm

Package contains:

- Single-use disposable SuperCross microcatheter
- Angled tip straightener
- Dispenser coil

Other items required but not provided:

- A guiding catheter with a minimum I.D. of 0.056" / 1.42mm fitted with a hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter $\leq 0.014"$ / 0.36mm
- 10ml syringe for flushing
- Sterile heparinized saline

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the SuperCross microcatheter packaging and components for damage.
2. Utilizing sterile technique, remove the SuperCross microcatheter from its packaging and transfer the dispenser coil to the sterile field.
3. Remove the catheter from the dispenser coil and wet the distal section with heparinized saline to activate the catheter's hydrophilic coating. Advance the angled tip straightener towards the distal tip such that it straightens the tip angle and maintains the tip in this straightened position in preparation for guidewire loading.
4. Inspect the SuperCross microcatheter for any bends or kinks. Thoroughly flush the SuperCross microcatheter with heparinized saline.

DEPLOYMENT PROCEDURE

Deploy the SuperCross microcatheter according to the following steps:

1. Backload the distal tip of the SuperCross microcatheter onto the proximal end of the exchange length guidewire that is already in place in the distal vasculature.
2. If the angled tip straightener has been used to straighten the tip of an angled SuperCross microcatheter tip to facilitate guidewire backloading, ensure that the guidewire is backloaded past the proximal edge of the angled tip straightener

3. If the angled tip straightener has been used and once the guidewire has been backloaded as described in step 2, carefully tear the angled tip straightener off the distal end of the catheter by gripping the proximal tab on the angled tip straightener and peeling it in the direction of the distal tip of the microcatheter. Discard the angled tip straightener once it has been removed from the catheter.
4. Under fluoroscopy, carefully advance the catheter until the marker at the distal tip is visible in the desired vascular space. **WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter or vessel perforation.**
5. If guidewire exchange is desired, after fluoroscopically confirming that the SuperCross microcatheter tip is positioned in the desired vascular space, slowly remove the guidewire.
6. Insert the desired exchange length guidewire into the luer-lock guidewire entry port. Advance the guidewire until it exits the catheter tip.
7. Under fluoroscopy, ensure that the guidewire is sufficiently advanced into the distal vasculature. Fix the position of the guidewire and slowly advance the SuperCross microcatheter into the desired vascular space.
8. If subselective infusion/deliver of diagnostic or therapeutic agents is desired, after fluoroscopically confirming that the SuperCross microcatheter tip is positioned in the desired vascular space, slowly remove the guidewire.
9. Create an air-free, fluid-to-fluid connection between the catheter hub and the injection syringe. Aspirate prior to injection, and observe the syringe to remove any air bubbles prior to injecting into the catheter.
10. Using the injection syringe, manually infuse the diagnostic or therapeutic agent. **WARNING: Do not exceed the maximum recommended injection pressure of 300 psi when injecting through the SuperCross catheter. Catheter damage or vessel injury could result.**
11. Ensure proper positioning of the exchange length guidewire prior to removing the SuperCross catheter and proceeding with the procedure.
12. Dispose of SuperCross microcatheter following standard hospital procedure.

STORAGE & HANDLING

No special storage or handling conditions.

REGULATORY INFORMATION

A pdf copy of this IFU is located at www.teleflex.com/IFU

This is the SuperCross Angled Tip microcatheter (Basic UDI-DI: **080190200000000000000291KX**) Summary of Safety and Clinical Performance (SSCP) location after the launch of the European Database on Medical Devices/EUDAMED: <http://ec.europa.eu/tools/eudamed>

For patient/user/third party in the European Union and in countries with identical regulatory regime (2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and go to your national authority. The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website: http://ec.europa.eu/growth/sectors/medical-device/contacts_en

LIMITED WARRANTY

Teleflex Medical LLC warrants that the SuperCross microcatheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Teleflex Medical LLC to be defective in workmanship or materials. Teleflex Medical LLC shall not be liable for any incidental, special, or consequential damages arising from the use of the SuperCross microcatheter. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of Teleflex Medical LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Teleflex Medical LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF TELEFLEX MEDICAL LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

SuperCross is a trademark of Teleflex Innovations S.à.r.l., or Teleflex Medical LLC, each a part of Teleflex Incorporated.

See the International Symbols Glossary on page 40.

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